Platelet PGD® Reagents

A. INTENDED USE

The Platelet PGD Reagents are only for use with the Platelet PGD Test.

B. SUMMARY AND EXPLANATION OF THE TEST

Platelet PGD Reagents are for use with the Platelet PGD Test. They are used in the processing of platelet samples prior to addition to the Platelet PGD Test Device for detection of aerobic and anaerobic Gram-positive and Gram-negative bacteria and in the processing of Platelet PGD Controls.

C. PRINCIPLES OF THE PROCEDURE

Refer to the Platelet PGD Test Instructions for Use (IFU).

D. REAGENTS AND MATERIALS

Materials Provided

Platelet PGD Reagents	REF P20R-CE	20 tests
	REF P100R-CE	100 tests

Includes the following:

	20 Test	100 Test
Reagent [1]	1 x 6 mL	2 x 12 mL
Reagent [2]	1 x 6 mL	2 x 12 mL
Reagent [3]	1 x 3 mL	2 x 6 mL

Reagents

Reagent [1]	Water, methanol and surfactants. Preservative: ProClin® 300
Reagent [2]	Water, sodium hydroxide and surfactants. Preservative: sodium azide
Reagent [3]	Tricine buffer with surfactants, anti-coagulants and protein (bovine, mouse, rabbit) stabilizers. Preservatives: ProClin 300 and sodium azide

See Reagent Precautions below.

Materials Required

Platelet PGD Test	REF P20CE	20 Tests
	REF P100CE	100 Tests

Materials Required But Not Provided - Refer to the Platelet PGD Test IFU.

E. WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use

Warnings

1. Read the Platelet PGD Test and Platelet PGD Controls IFU completely before using the product. Follow the instructions carefully. Not doing so may result in inaccurate test results.
2. The Platelet PGD Reagents have been validated only for use with the Platelet PGD Test and Platelet PGD Controls.
3. Each operator performing the test must be able to distinguish between the following colors: Green, Yellow, Blue and Red.
4. Do not use materials after their stated expiration dates.

Reagent Precautions

Reagents were classified according to applicable GHS and European Community (EC) Directives. Applicable Classification, Hazard (H) and Precautionary (P) statements are listed below. Refer to SDS for complete Precautionary Statements.

Reagent [1] contains methanol and the signal word is Danger.

H226 Flammable liquid and vapor.

H302 Harmful if swallowed.

H312 Harmful in contact with skin

H319 Causes serious eye irritation

H370 Causes damage to organs.

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P312 Call a Poison Center or doctor if you feel unwell.

P501 Dispose of contents/container in accordance with local and national regulations.

Reagent [2] contains sodium hydroxide and the signal word is Warning.

H315 Causes skin irritation

H319 Causes serious eye irritation

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351 If in EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical attention.

P501 Dispose of contents/container in accordance with local and national regulations.

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Reagent [3] contains n-Dodecyl-N,N-dimethyl-3-ammonio-1-propanesulfonate (DDAPS) and the signal word is Warning.

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P305+P351+ If in EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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Reagents 2 and 3 contain sodium azide. Contact with acids liberates very toxic gas.

General Safety Precautions

Follow good laboratory practices and use Universal Precautions when handling all samples and materials.^{1,2} Dispose of all test materials as bio-hazardous waste according to your laboratory procedure and applicable regulations.

Handling Precautions

Handle and perform test properly:

1. Do not combine leftover volumes of Reagents 1, 2, or 3.
2. Do not remove dropper tips from bottles.
3. Do not touch exposed dropper tips.
4. Recap bottles immediately after use. Do not interchange bottle caps. The cap color must match the label color.
5. Do not use test components beyond the expiration dates printed on the labels. Always check expiration dates prior to performing test.
6. Do not use Reagents if they have not been properly stored at 2 - 8 °C. It is not necessary to equilibrate Reagents to room temperature prior to use.

F. STORAGE INSTRUCTIONS

1. Store Platelet PGD Reagents at 2 - 8 °C. Once opened, use prior to the expiration date on the v

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