Platelet PGD® Controls

A. INTENDED USE

Platelet PGD Controls are only for use with the Platelet PGD Test as assay Quality Control Samples to verify the performance of the Platelet PGD Test.

B. SUMMARY AND EXPLANATION

The Platelet PGD Controls are used to verify the user's ability to properly perform the test and interpret test results. Failure of the Platelet PGD Controls to perform as expected may also indicate deterioration of the Platelet PGD Reagents or Platelet PGD Test Device.

C. PRINCIPLE OF THE PROCEDURE


The Positive Control contains bacterial antigens (in a matrix containing a buffered solution of animal protein, human platelet lysates and preservatives) that react with all four (4) detection lines on the Platelet PGD Test Device. The Negative Control contains the same matrix, without bacterial antigens, and therefore will not react with any of the detection lines.

Two drops of Positive or Negative Control are added to a Microfuge Tube then added with the Platelet PGD Test. Reagents 2 and 3 are then added, the sample is mixed, and the entire volume is added to the Platelet PGD Test Device. After the test has run to completion, the results are analyzed. The Positive Control should result in four visible detection lines, two in the Gram-positive (GP) and two in the Gram-negative (GN) Test Result Windows. The Negative Control should result in no visible detection lines.

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

CONTROL - Phosphate buffered saline, platelet lysates and protein (human, rabbit) stabilizers. Preservatives: ProClin 300® and sodium azide

CONTROL + Phosphate buffered saline, Lipoteichoic acid, bacterial antigens, platelet lysates and protein (human, rabbit) stabilizers. Preservatives: ProClin 300 and sodium azide

Materials Required

Platelet PGD Test	REF P20CE	20 Tests
	REF P100CE	100 Tests

Materials Required But Not Provided - Refer to the Platelet PGD Test Instructions for Use (IFU).

E. WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use

Warnings

-  The Platelet PGD Controls contain human sourced and/or potentially infectious components. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents be handled in accordance with the OSHA Standard on Bloodborne Pathogens using Universal Precautions.¹ Bio-safety level 2 or other appropriate bio-safety practices should be used for materials that contain or are suspected of containing infectious agents.²
- The human derived components within the Positive and Negative Controls are non-reactive for HBsAg, HIV-1 RNA, anti-HIV-1/HIV-2, anti-HCV, HCV RNA, anti-HTLV-I/II and West Nile Virus (WNV) RNA.
- Read the Platelet PGD Test IFU completely before using the product. Follow the instructions carefully. Not doing so may result in inaccurate test results.
- Do not use materials after their stated expiration date.

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Reagent Precautions

- This product contains sodium azide as a preservative. Contact with acids liberates very toxic gas. For product not classified as dangerous per European directive 1999/45/EC as amended, Safety data sheet available on request.
- This material and its container must be disposed of as bio-hazardous waste according to your laboratory procedure and applicable regulations.

Handling Precautions

- Handle Controls properly:
- Do not combine partially used bottles of Positive or Negative Control.
 - Do not remove dropper tips from bottles.
 - Do not touch exposed dropper tips.
 - Recap bottles immediately after use. Do not interchange bottle caps. The cap colors must match the label colors.
 - Check expiration dates prior to performing test.
 - Do not use if Controls have not been stored at 2 - 8 °C. It is not necessary to equilibrate Controls to room temperature prior to use.

F. STORAGE INSTRUCTIONS

Store Platelet PGD Controls at 2 - 8 °C. Once opened, use prior to the expiration date on the vial.

G. INDICATIONS OF INSTABILITY

- Inspect Control bottles for precipitate. Do NOT use if precipitate is present.
- Failure of the Platelet PGD Controls to perform as expected may indicate deterioration of the Platelet PGD Reagents or the Platelet PGD Test Device.

H. CONTROL TEST PROCEDURE

Refer to the Platelet PGD Test IFU for *Control Processing* and *Performing the Test*.

I. INTERPRETATION OF RESULTS

Refer to the Platelet PGD Test IFU for *Interpretation of Results*.

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K. REFERENCES

- 29CFR 1910.1030; Occupational Safety and Health Standards: Bloodborne Pathogens
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization; 2004

Verax Biomedical and PGD are trademarks of Verax Biomedical Incorporated.

ProClin 300 is a trademark of Rohm and Haas Company.

Keys to Symbols Used

IVD	<i>In vitro</i> diagnostic medical device
LOT	Lot Number
REF	List Number
	Expiration Date
CONTROL -	Negative Control
CONTROL +	Positive Control
	Temperature limitation
	Manufacturer
EC/REP	Authorized Representative
	Attention, see instructions for use

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